

# RB-208/210 OPERATOR'S MANUAL

Published: 7-22-03

Revision: 1.0

Phone: (216)292-7373

E-mail: rexonsales@rexon.com Visit us at www.rexon.com

Fax: (216)292-7714



# MODEL RB-208/210 Series Integrated PMT Base

The Rexon Model RB-200 series Integrated PMT Base provides a simple solution for a reliable interface to photomultiplier tubes (PMT). The unit provides a high speed, low noise preamplifier, regulated, low ripple high voltage power supply and divider network is suitable for use in a wide variety of scintillation detector applications. The low noise design of the HV supply eliminates spurious counts often associated with lower quality units while eliminating the need for a costly external supply and associated cabling The PMT socket is a standard 14 pin configuration and is available for both 8-stage and 10-stage dynode configurations.

The RB-200 features a convenient, sub-miniature LCD that displays the HV setting. The HV supply can be adjusted either the locally or remotely. In addition, a logic level output signal indicates HV normal or fault status. The RB-200 features high stability and low ripple PMT bias supply necessary to meet the most demanding requirements.

The Integrated PMT Base provides an interface for use with a variety of MCA cards allowing remote control of the HV, monitoring of the HV supply status and low voltage supply direct from the MCA card.

Since the RB-200 provides all of the functionality necessary to properly bias the PMT, condition and drive the signal system integration becomes easy. Connection to the unit requires simply a low voltage power converter and a coax cable with a BNC connector. An internal jumper allows output impedance to be set for 50 ohms or 93 ohms to match the particular coax cable used.

The RB-200 is available in several versions. It can be supplied in either a 8-stage dynode or 10-stage dynode configuration. A preamplifier only version provides an SHV connector when external HV is used.

Both 120VAC/60Hz and 220VAC/50Hz power converters are available as well as a wide variety of international converters.

Phone: (216)292-7373

E-mail: rexonsales@rexon.com Visit us at www.rexon.com

Fax: (216)292-7714

Consult our factory for other custom versions, configurations or gain selections.





Figure 1: RB-200 SERIES CONNECTOR PANEL

# **Specifications**

Display: 3.5 digit LCD

Display resolution: 1V

Voltage range: 0 to 2000VDC (positive)

Operator controls: HV adjustment, remote/local selection

(internal)

Signal connector: BNC type
Output impedance: 50 or 93

Output impedance: 50 or 93

HV output current: 0 to 0.5mA

HV ripple (, full load, P-P): <0.0005%

HV stability: <0.005%/hr

HV Temperature coefficient (typ.): <25ppm/°C

Operating temperature: -10° to +50°C (85% RH max.)

Storage temperature: -25° to +95°C (95% RH max.)

Power requirements: 120VAC/60Hz to 12VDC @ 300mA power converter (220VAC/50Hz

available)

Rexon Components, Inc. 24500 Highpoint Road Beachwood, Ohio 44122-1022 Phone: (216)292-7373 Fax: (216)292-7714 E-mail: rexonsales@rexon.com Visit us at www.rexon.com



Direct from MCA (+12V)

Phone: (216)292-7373

E-mail: rexonsales@rexon.com Visit us at www.rexon.com

Fax: (216)292-7714

Dimensions (Dia. x H): 2.2" x 5.4" (5.6 cm x 13.7 cm)

Weight: 0.8 lb (0.28 kg)

# Set-up and connections:

All user connections are located on the RB-200 connector panel. The following section describes these connections.

## POWER:

The RB-200 requires a 12VDC @ 300mA power supply which is supplied with the unit. The standard DC power jack (positive center pin) plugs directly into the base of the unit. There is no power switch, therefore the RB-200 turns on immediately when the power cord is connected.

## SIGNAL OUTPUT:

The output signal from the RB-200 is available on the BNC connector. This signal is 50 ohm impedance and suitable for direct connection to most Multi-Channel Analyzers (MCA) or counters and a variety of other data acquisition systems. This output signal is positive polarity and its amplitude is controlled by adjusting the PMT high voltage supply. Standard RG-58 coaxial cable is normally used for this connection.

## **HV ADJUSTMENT:**


The potentiometer on the connector panel allows the HV to be set by the user. The HV setting is indicated on the display (LCD).

CAUTION: It should be noted that some PMTs high voltage ratings that are considerably lower than the 2000V range of the HV supply. It is necessary that the user be aware of the maximum rating for the PMT used and that when making adjustment, be cautious not to exceed that rating.

### MCA CONNECTOR:

The MCA connector is provided for applications where it is desirable to supply power and input connections directly from the MCA. Although some MCA's provide this capability, others do not. Refer to the appropriate MCA Operator's Manual for interface requirements. As noted, the RB-200 requires 12VDC @ 300ma.





#### MCA SIGNAL DESCRIPTION:

## **HV CONTROL:**

This input signal, a positive, 0-5V signal, applied externally, via the MCA, to adjust the high voltage supply, thereby altering the overall gain of the unit. Internally, the REM/LOCAL switch must be placed in the REM (remote) position. This 0-5V signal controls the high voltage supply over it's entire, 0 to 2000V range.

CAUTION: It should be noted that some PMTs high voltage ratings that are considerably lower than the 2000V range of the HV supply. It is necessary that the user be aware of the maximum rating for the PMT used and that when making adjustment, be cautious not to exceed that rating.

#### **HV STATUS:**

This signal is a digital (0-5V) output indicating whether the HV supply greater than approximately 200VDC. This signal is pulled to ground when the HV is greater than 200V and is pulled high (10K ohms /+5V) when the supply has failed (below 200 VDC).

# Operation:

After the appropriate connections have been made, plug in the power cord and apply power to the unit. The HV adjustment may be used to set the gain of the system. The user can verify that the gain is appropriate and that the signal amplitude is sufficient by observing the signal output on an oscilloscope. A combination of HV adjustment, to obtain optimum performance from the PMT and external gain adjustment, such as MCA input amplification, will provide a suitable spectrum for analysis purposes. In counter applications, the output signal needs to be high enough to exceed a threshold adjustment generally provided by the counter. Note, the output signal is an analog signal and as such, depending on gain and source energies, can vary greatly in amplitude. **Do not input directly to a device that requires a digital signal.** 

Phone: (216)292-7373

E-mail: rexonsales@rexon.com Visit us at www.rexon.com


Fax: (216)292-7714

REXON
Components, Inc.

24500 Highpoint Rd. Beachwood, Oh 44122 Phone: (216) 292-7373 Fex: (216) 292-7714

Job No.	9671
Customer	Anjali Irc.
Type 'PV'	TM51 PX 10000 /2-883
Date	1-28-03
Tested By	KE
Serial No	030/28-1

Aualyser Potce Spile	
Conv Gain 1024	
DOAMP Sals MCA Card Bigolas 10x1 (RB210)	YES
PMT REXOD 800P (WE7878)	Z DNO
Source CS137 Sold applied 6 from PMT	<b>Y</b>
With Voltage 12 CC	103




REXON
Components, Inc.

24500 Highpoint Rd. Beachwood, Oh 44122 Phona: (216) 292-7373 Fax: (216) 292-7714

Job No.	91071
Custemer	Anjali Inc
	M51 PX1000m / 28821
Date	1-28-03
Tested By	KE
Serial No	030/28-1

	Analyser Apter sollo
	Conv Gain 1024
,	Amp Sain McAcard & polar 10x4 ( B. 8210)
	PMT Rexon 800P (WE7878)
	Source CS 137 506" Appen 18" from PINT
	Ulah Valtona 12.00





Components, Inc. 24500 Highpoint Rd. Beachwood, Oh 44122 Phone: (216) 292-7373 Fax: (216) 292-7714

Aptec 5011 Analyser Job No. Conv Gain 024 Customer AB210 Type PYT m51 PX 1000m 12 PRIMATIO GAIN MCA Cord PINT Date **Tested By** High Voltage 1200 Serial No 030/28

☐ YES

□ NO

